

How to Produce and Use the siRNA Cocktails for RNAi Knockdown

QUAN KANG 11/3/2003; Edited by T.-C. He

dsRNA Transcription In Vitro Using T7 RiboMAX™ Express Large Scale RNA Production System

1. Prepare the linear DNA template.
2. Assemble transcription reaction at room temperature.

T7 Reaction Components	Sample Reaction
RiboMAX™ Express T7 2X Buffer	10.0μl
Linear DNA template(0.5-1μg total)	1-8μl
Nuclease-Free Water	0-7μl
Enzyme Mix, T7 Express	<u>2.0μl</u>
Total Volume	20μl

3. Mix gently and incubate at 37°C (water bath) for 1-2 hrs.
4. Incubate at 75°C for 5 minutes then leave the mixture on the hood to cool to room temperature.
5. Add RQ1 RNase-Free DNase to a concentration of 1unit(1unit=1μl) per microgram of template DNA to remove the DNA template and unincorporated rNTPs .
6. Incubate for 15 minutes at 37°C.
7. Extract dsRNA by PC-8.
8. Precipitate dsRNA by 100% ethanol and wash it by 70% ethanol twice.
9. Dissolve the dsRNA in nuclease-free water to a volume identical to that of the transcription reaction, store at -80°C.

dsRNA Digestion Using ShortCut RNase III (NEB) and siRNA Cocktail Purification

1. Transcribe and purify dsRNA in vitro.
2. Assemble RNase III digestion reaction in a 1.5ml microfuge tube.

ddH ₂ O	75-x μl
10×ShortCut Reaction Buffer	10μl
dsRNA	x μl
ShortCut RNase(NEB)	5μl
10×Mncl ₂	<u>10μl</u>
total volume	100μl

3. Mix and incubate at 37°C (water bath) for 20 minutes.
4. Add 10µl 10×EDTA to stop the reaction.
5. Remove the product to a Microcon G-30 tube, centrifuge at 14,000g for 15 minutes. Collect the siRNA in the tube, store at -80°C.

Tansfect siRNA Into Mammalia Cell Lines with Lipofectamine

For 24-well plate

1. Plate the cells in 24-well plate at 8-12hrs prior to transfection (approx 30-50% confluent).
2. prepare the Lipofectamine/siRNA mixture(Per well):
100µl serum free DMEM + 3µl Lipo + 5~20µl siRNA
3. Mix well and incubate at room temperature for 20-30 minutes.
4. Meanwhile, remove the medium from every well, wash cell gently with 0.5-1ml of serum free DMEM.
5. Add 200µl serum free DMEM to each well, and return plate to incubator for 5-10 minutes.
6. Add the Lipo/siRNA mixture to each well, rock gently then return it to incubator.
7. After 3-5hrs, remove the medium and replace with 1.0ml complete medium.

Note: the amount of siRNA is up to your experimental purpose.

For silencing the endogenous genes without any stimulation, try 5~10µl siRNA, but for silencing genes after stimulation(eg. After infection), try 10~20µl.