

Mouse Tail DNA Extraction for Genotyping

(Updated by Guo-Wei Zuo, 09/26/2009)

1. Label 1.5ml tubes with appropriate animal numbers.
2. Cut a small piece (approx. 1 mm³) of mouse tail, using a scalpel. Samples can be stored at -70°C if needed.
3. Add 100ul Proteinase K buffer, vortex briefly, and incubate overnight at 50°C.

Prot K digestion buffer:

10mM Tris pH 8
100mM NaCl
10mM EDTA
0.1mg/ml Prot K
0.5% SDS

4. Add 200ul ddH₂O, Cfg for 1min at top speed to bring down moisture on sides of tubes.
5. Pour supernatant to a new tube, add 250ul PC-8 (phenol:chloroform), vortex, cfg for 2min at top speed.
6. Transfer the aqueous phase to a clean/labelled 1.5ml tube.
7. Add 700ul cold ethanol and invert to mix. (DNA ppt may not be visible.)
8. Cfg in microfuge for 5min at RT.
9. Aspirate supernatant; add 500ul of 70% ethanol, vortex well, and spin down for 1 min.
10. Aspirate supernatant; Spin down for 10 sec and aspirate the residual liquid completely. Pellets may not be clear.
11. Add 200ul ddH₂O to each pellet and vortex briefly.
12. Store samples at -20°C.