


Extract DNA from Agarose Gel Using Montage DNA Gel Extraction Device (Millipore)

Quan Kang 11/5/2003; Updated by JK Park 4/20/04, Commented by TCH

1. Prepare a **six-well** minigel of **0.5% to 0.8%** agarose [**Note: The gel is thicker than regular ones, and you will need 70-100ml to pour such a gel**].
2. Load the DNA digestion mix (usually **100ul** reaction plus **20ul** of 6x GSB dye) [**Note: If you have spare wells, leave at least one blank well between samples in order to minimize potential cross-contaminations**].
3. Run the gel at 60 to 75 voltages for 30 min. to 60 min [**Note: The running time depends on the size of your interested DNA fragment**].
4. Locate the DNA band of interest in the gel using a UV lamp and cut the gel containing the DNA. Trim down the gel piece as much as possible. **The gel slice should be under 100 μ l in volume or 100mg in mass.**
5. Make sure the montage DAN gel extraction device (a.k.a. ULTRAFREE-DA spin column, Millipore) is assembled as follows:



6. Place the gel slice in the Gel Nebulizer and seal the device with the cap attached to the vial.
7. Spin the assembled device for 10-15 minutes at 5,000 \times g (~7,500 rpm).

8. Discard the Sample Filter Cup and Gel Nebulizer units.
9. Recover the eluted DNA, and bring the volume to 200ul with ddH₂O. Add 100ul 7.5M ammonium acetate.
10. Extract the DNA mix with phenol-chloroform by adding 250ul PC-8, vortexing, and briefly spinning down [**Note: Repeat the PC-8 extractions if dirties still appear at the interphase**]. Transfer the clear DNA solution to a new 1.7ml microfuge tube.
11. Precipitate DNA by adding 4ul seeDNA and 600ul to 700ul of ice-cold 100% ethanol. Centrifuge the sample at top speed for 5 min. Wash the pellet with 70% ethanol twice.
12. Air-dry the pellet. Dissolve the DNA in **12ul** ddH₂O, which is ready for use (e.g., use **3ul** per ligation reaction).