

How to Use Bio-Rad Semi-Dry Transfer Apparatus for Western Blotting Analysis

(Quan Kang 8/28/03, Commented by TCH)

1. Run SDS-PAGE pre-cast gels (e.g., from Bio-Rad or ISC BioEXpress) as usual.
2. Prepare transfer membrane: cut the PVDF membrane to the dimension of the gel. Wet the membrane by soak it into 100% methanol 10~15 sec, then wash it with distilled water 3 times, keep the membrane in Transfer Buffer (**Note: You must first wet the membrane in 100% methanol and never let it dry afterward**).
3. Remove safety cover and stainless steel cathode. Assemble the transfer sandwich in the following order (from the bottom platinum anode up): one piece of pre-wet extra thick filter paper (or 3~6 general thick filter papers) >>> Pre-wet PVDF membrane >>> SDS-PAGE gel (**Optional: the gel can be pre-equilibrated in the Transfer Buffer for 5-20min**) >>> one piece of pre-wet extra thick filter paper (or 3~6 general thick filter papers). **Note: Roll out air bubbles between layers with a 5ml pipette after assembly the sandwich.**
4. Secure top stainless steel cathode and safety cover. Start the transfer using the following suggested conditions: **Mini gel**, 10-15V, 250mA for 20min (high MW proteins, run 25~30min); **Large gel**, 25V, 250mA for 30min.
5. Turn off the power supply, unplug electrodes, and remove blot. Complete Western blotting by following the conventional Western Blotting protocol.

Tips:

- 1) The PVDF membrane and filter papers must be 100% equilibrated.
- 2) The filter papers should be cut to the exact size of the gel (this forces the current to flow only through the gel and not through overlapping filter paper).
- 3) Any bubbles between the layers will block current flow and prevent protein transfer.
- 4) Once assembled, do not move the top stainless cathode, any shifting of the transfer stack after assembly will distort the transfer pattern.

- 5) Multiple gels can be transferred, simply put a sheet of porous cellophane or dialysis membrane equilibrated with Transfer Buffer between each transfer sandwich.

Transfer Buffer (for PVDF membrane)

800 ml Methanol (final concentration = 20%)

12.12g Tris Base

57.63g Glycine

Add ddH₂O to 4 Liters

(Note: pH should be 8.3~8.6, but NEVER add acid or base to adjust it).