

## MOlab Standard Protocol for Preparing RT-PCR cDNA Stocks for qPCR

(Lan Zhou 11/02/00, updated by Yang Bi 02/24/08, commented by TCH)

### NOTE:

- A) Yang has recently compared the RT enzyme from NEB (**Cat# M0253S**) vs. Invitrogen's Superscript II, and found no significant difference between those enzymes. Thus, the **NEB RT enzyme is the recommended choice** in our lab. **However, we do not use any manufacturer's protocol!**
- B) After numerous pre-testing runs by former/current lab members, we decide that the RT enzyme in each reaction can be reduced to as low as **0.25 $\mu$ l** (the original suggested volume was 1.25 $\mu$ l).
- C) These conditions work well for most regular RT-PCR cDNAs prepared for real-time PCR assays, especially if your RNA is isolated from cultured cells and/or fresh tissue samples. However, if you are concerned about the yield and quality of your RNA samples, you are always recommended to check a couple microliters of your RNA samples on RNA gel (see **Protocol C12**).
- D) Make sure you use enough total RNA for each RT reaction. Typically, a subconfluent T-25 flask yields about 5-10 $\mu$ g total RNA, and a subconfluent **100mm dish yields about 10-20 $\mu$ g total RNA**. Please keep in mind that the RNA yield will be much lower if the cells are cultured in low serum medium 8-12h prior to RNA isolation.
- E) In our TRIzol or Nucleozol RNA isolation protocol, we do not routinely use RNase-free DNase to digest away potential genomic DNA. Thus **-RT tube** is set up in some cases. However, we have found that **the -RT tube is not necessary in most cases**; so **it is optional**.

1. Turn on two heating blocks: one at 70°C and one at 37°C.

2. Make Hexamer mix [**prepare half of the volume if -RT tube is not set up**]:

	1X
Hexamer(random primer, 1 $\mu$ g/ $\mu$ l)	4 $\mu$ l
SSB (0.5 $\mu$ g/ $\mu$ l) ( <b>optional</b> )	2 $\mu$ l
RNase-free H2O	3 $\mu$ l
<b>Total volume</b>	9 $\mu$ l ( <b>or 4.5<math>\mu</math>l per RT</b> )

3. Make Hexamer mix-RNA [**prepare half of the volume if -RT tube is not set up**]:

Hexamer mix	9 $\mu$ l
Total RNA (10 $\mu$ g) or mRNA (0.5 $\mu$ g)	x $\mu$ l
RNase-free H2O	xx $\mu$ l
<b>Total volume</b>	25 $\mu$ l ( <b>or 12.5<math>\mu</math>l per RT</b> )

4. Incubate @ 70°C for 3-5 minutes.

5. Prepare RT mix [**prepare half of the volume if -RT tube is not set up**]:

5 X First strand buffer	10 $\mu$ l
0.1M DTT	4 $\mu$ l
10mM dNTPs	2 $\mu$ l
RNasin ( <b>optional</b> )	0.4 $\mu$ l
<b>Total volume</b>	16.4 $\mu$ l ( <b>or ~8<math>\mu</math>l per RT</b> )

6. Prepare +/- RT Reactions:

	<b>+RT tube</b>	<b>-RT tube (optional)</b>
Hexamer mix-RNA	12.0 $\mu$ l	12.0 $\mu$ l
RT mix	7.8 $\mu$ l	7.8 $\mu$ l
<b>RT enzyme (NEB)</b>	0.25 $\mu$ l	0.0 $\mu$ l
RNase-free H2O	0.0 $\mu$ l	0.25 $\mu$ l
<b>Total volume</b>	<b>20<math>\mu</math>l</b>	20 $\mu$ l

7. The cycling program of RT-PCR:

37°C X 60 minutes

95°C X 1 minute (for killing RT; RT interferes with Taq DNA Pol)

4 °C Hold or kept at -80°C freezers

8. Add 80µl ddH<sub>2</sub>O to get 100µl original cDNA, and then take 10µl out, add 40µl ddH<sub>2</sub>O to dilute into 5 times (**total volume = 100µl, aka., RT-PCR cDNA stock**), and then use the diluted cDNA for real-time PCR. **It's a good practice to aliquot the initial 100µl into multiple aliquots. Keep the aliquots at -80°C.**
9. For most real-time PCR reactions, the cDNA mix needed to be further diluted 5 to 100 times, depending on transcript abundance of the gene of your interest. **The Ct or Cq values for internal reference genes should be ideally between 15~20 cycles.**