

T4 DNA POLYMERASE BLUNT-ENDING/FILL-IN REACTION FOR CLONING

Updated by Yi Shu 05/21/2017; edited by TCH

- 1) The restriction enzyme digested DNA should be precipitated with ethanol and washed with 70% ethanol twice;
- 2) Set up a 100µl reaction as the follows:

10 x NEB Buffer 2.1(NEB)	10.0
dNTPs (10mM each)	10.0
ddH ₂ O	79.0
Digested pDNA	dried
T4 DNA Polymerase (NEB)	1.0
	<hr/>
	100.0 µl

- 3) **Incubate at 37°C for 15min ~ 30min.**
- 4) Precipitate blunt-end DNA with ethanol (see Ethanol Precipitation protocol) or directly load the reaction mix to 0.5% gel for purification.

- NOTE:**
- 1) T4 DNA polymerase exhibits both 3'- and 5'exonuclease activities so it can fill in or chew back either 3'- or 5'-protruding ends.
 - 2) To avoid unnecessary chewing back, one must ensure a large excess of dNTPs in the reaction.
 - 3) When two different restriction sites are blunt-ended and ligated, there is a 25% chance restoring either of the restriction sites.
 - 4) For fill-in reactions only: T4 DNA Polymerase can be used in NEBuffers 1.1, 2.1, 3.1 and CutSmart® Buffer as well as NEBuffers 1-4 and T4 DNA Ligase Reaction Buffer.
 - 5) For blunting reactions requiring removal of overhangs: T4 DNA Polymerase can be used in NEBuffers 1.1, 2.1, and CutSmart Buffer as well as NEBuffers 1, 2, and 4 and T4 DNA Ligase Reaction Buffer. NEBuffers 3.1 and 3 are not recommended when overhang removal is required.
 - 6) **Elevated temperatures, excessive amounts of enzyme, failure to supplement with dNTPs or long reaction times will result in**

recessed ends due to the 3' → 5' exonuclease activity of the enzyme.