

ThermoSequenase DNA Cycle Sequencing Protocol  
Adapted from BV's Cookbook by T.-C. He, 10/24/2000

This protocol describes the use of the Amersham DyeTerminator/ThermoSequenase Kit for Cycle sequencing PCR products, plasmids, or BACs. See [http://www.apbiotech.com/technical/documentation/dna\\_sequencing/XY0430003.pdf](http://www.apbiotech.com/technical/documentation/dna_sequencing/XY0430003.pdf) for more information.

**A. Prepare templates:**

**A-1. Plasmid DNA minipreps:**

Plasmid DNA purified 2ml overnight coli culture using alkaline lysis procedure can be routinely used in cycle sequencing reactions. If the final volume of a miniprep is 50-70ul, it usually requires 4ul per primer reaction (or 1ul per A/C/G/T reaction).

**A-2. BAC clones:**

If BAC DNA is used, it usually requires approximately 0.5-1.0ug per sequencing reaction.

**A-3. PCR products:**

If you are sequencing a PCR product, you may want to remove excess amplification primers from the PCR reaction. Purify the PCR product by either of the methods listed below:

**I) NaClO<sub>4</sub> Precipitation:**

- 1). Bring volume of PCR reaction to 600 ul in TE or ddH<sub>2</sub>O.
- 2). Perform PC8 extraction.
- 3). Transfer supernatant to a new 1.7ml tube, add:

2 ul	glycogen (or 3ul seeDNA)
200 ul	2M NaClO <sub>4</sub>
400 ul	2-propanol
- 4). Microfuge 10', RT.
- 5). Decant supernatant, wash twice with 70% EtOH.
- 6). Air dry the pellet, and resuspend DNA pellet in appropriate volume of ddH<sub>2</sub>O to allow convenient delivery of about 0.5-5ug PCR product per sequencing reaction (usually about 20ul total volume).

**II) Gel Purification:**

- 1). Run PCR product on agarose gel of appropriate concentration (0.7-2.0%).
- 2). Cut out appropriate band.
- 3). Place gel slice into upper reservoir of Spin-X tube (Costar), or regular spin-column.
- 4). Microfuge 15', RT.
- 5). Bring volume of liquid collected in lower reservoir of Spin-X tube to 400ul.
- 6). Add:
 

2 ul	3ul see DNA
40 ul	7.5M NaOAc
800 ul	100% EtOH
- 7). Microfuge 5', RT.
- 8). Decant supernatant, wash twice with 70% EtOH.
- 9). Air dry and resuspend DNA pellet in appropriate volume of TE or ddH<sub>2</sub>O to allow convenient delivery of about 0.5ug PCR product per sequencing reaction (usually about 20ul total volume).

## **B. Cycle Sequencing Reactions**

1. Set up Reaction Mixture:
 

2ul	10 X Reaction Buffer
4ul	plasmid DNA, minipreps
	[or ~2 ul PCR product (50-500 ng) ]
	[or for BAC DNA, use ~10% of a Nucleobond Midiprep]
1ul	Primer (20-40 ng/ul; or 0.5-2.5 pmol)
12ul	Water (or 20 ul minus the volume of all other components)
1ul	ThermoSequenase (USB, 4U/ul) (in some cases, 2ul can be used)
2. Make 4 master termination mixtures (G, A, T, C) containing 2N ul dGTP termination mix and 0.5N ul of the appropriate [ $\alpha$ -33P]ddNTP, where N = the total number of sequencing reactions to be run. (Alternatively, use dITP termination mix instead to dGTP termination mix; dITP is recommended for resolving compressions, while dGTP supposedly yields more even band intensities).
3. Aliquot 2.5ul of each master termination mixture (G, A, T, C) to separate wells of a 96-well PCR plate.

4. To each well, add 4.5 ul of Reaction Mixture. Overlay with 10-20ul mineral oil.

6. Cycle reactions:

<b>dGTP</b>	<b>Reactions</b>	<b>dITP Reactions</b>
Stage 1	94°C × 2' ×1 cycle	Stage 1 94°C × 2' ×1 cycle
Stage 2	92°C × 20"	Stage 2 92°C × 20"
	55°C × 30" × 25 cycles	47°C (orTm-8)× 30" × 25cycles
	72°C × 1'	60°C × 4'

25 cycles has been successful with fairly robust reaction products. More cycles may increase background. Note:  $T_m = 4(G+C) + 2(A+T)$ .

7. Add 4 ul **Stop Solution** to each reaction. Check for blue dye in the aqueous layer.

8. Heat to 70°C for 5'. Load 3-4 ul, avoiding oil layer.

**NOTE:** DMSO can enhance sequencing with high GC templates. The protocol is the same except that the following Reaction Mixture is substituted for that in step 1 above:

2 ul	10 X Reaction Buffer
4 ul	Miniprep plasmid DNA
1 ul	Primer (20-40ng/ul; 0.5-2.5 pmol)
11.5 ul	ddH <sub>2</sub> O
0.5 ul	DMSO
1 ul	ThermoSequenase (USB, 4U/ul)