

Alizarin Red S Staining of Cells Cultured in Mineralization Medium

(By Hongwei Cheng @ 06/11/01; updated by Yi Zhu @ 01/08/2024; Commented by TCH)

A. Chemicals/Reagents for Mineralization Medium

- 1). L-Ascorbic acid (AA) ($C_6H_7O_6Na$)
- 2). Sodium β -glycerophosphate pentahydrate (bGP)($C_3H_{19}Na_2O_{11}P$), MW 308.13 g/mol

B. Stock Preparation

- 1). **L-Ascorbic acid (AA) @ 0.25g/ml (= 500x stock)**
1.25g Ascorbic acid in 50ml 1x PBS; syringe-filtered and aliquoted/kept @ -20°C.
- 2). **Sodium β -glycerophosphate pentahydrate (bGP) @ 1.0M (= 100x stock)**
3.081g β -glycerophosphate in 10ml 1x PBS; syringe-filtered and aliquoted/kept @ -20°C.
***NOTE:** It is more desirable to mix AA and bGP solutions before use.*

C. Working/Cell Culture Medium Preparation

- 1). The total volume varies depending on experimental design.
- 2). For example, to prepare 50ml AA/bGP-containing working medium, add 100ul of **AA (500x stock)** and 500ul of **bGP (100x stock)** to 49.4ml **Complete DMEM** (in a 50-ml conical centrifuge tube) [*The final concentration for ascorbic acid is 50 μ g/mL, and β -glycerophosphate 10mM, respectively*].

D. Cell Culture in Mineralization Medium

- 1). Seed cells in cell culture plates (24-, 12, or 6-well) for 4-6h to allow full attachment.
- 2). Replace the medium with the AA/bGP-containing mineralization medium prepared above.
- 3). You can add other treatment at this stage as well (e.g., adenoviral infection, etc.).
- 4). Depending on cell types and/or experimental design, keep the culture in a CO₂ incubator for 5-14 days, and proceed to Alizarin Red S staining (*see next section*).

***(NOTE:** The culture medium may turn to yellowish color after 3-5 days. This is normal and very common. You **should NOT ever change the medium**. Instead, you can **add some fresh complete DMEM medium**, and let the culture continue till the endpoints of the assay).*

E. Alizarin Red S Staining

- 1). Prepare **2% Alizarin Red S Solution** (usually stable for at least 2-3 weeks): Dissolve 2.0g of Alizarin Red S (Sigma-Aldrich) in 100ml ddH₂O. Mix well (*using magnetic stirring bar if necessary*)[***NOTE:** Use Whatman filter paper to remove debris or precipitates if necessary*]
- 2). Prepare **1x PBS pH 4.2 solution**: To 100ml 1x PBS, add 10% (v/v) HCl dropwise to adjust pH to 4.2 (*approx. 500 μ l of 10% HCl for 100ml 1x PBS*).
- 3). Fix cells in 2.5% glutaraldehyde (freshly prepared in 1x PBS buffer) for 10-15min at room temperature.
- 4). Remove fixation solution and wash the fixed cells with 1x PBS pH4.2 twice.

- 5). Add filtered 2% Alizarin Red Solution to the fixed cells and leave in 37°C incubator for 10-20min (Note: monitor staining under microscope every 2-5 min.).
- 6). Remove Alizarin Red Solution and rinse the cells with regular 1x PBS 1-2 times.
- 7). Record staining results with a bright field microscope (use non-phase contrast; positive stained nodules are in orange red).

GENERAL COMMENTS

- 1). It is important to wash the fixed cells with PBS adjusted pH to 4.2 to get brilliant red color.
- 2). **Alizarin Red can be quantified by dissolving the stain in a solution of 20% methanol and 10% acetic acid in water for 15min, followed by OD reading on a spectrophotometer at 450nm.**
- 3). Staining should be stable for 1-2 weeks in either dry or wet condition, but the best images are obtained immediately after staining.

SAMPLE STAINING RESULTS

