

# NucleoZOL Total RNA Isolation Protocol

(Xinyi Yu & TCH @ 08/21/2017)

**NucleoZOL Reagent, Cat# 740404.200** (BD ClonTech-Takara)

- 1) **LYSE CELLS:** For cells grown in monolayer, remove cell culture medium and lyse cells by addition of **1mL of NucleoZOL** to the culture dish (e.g., **1ml per one confluent 100mm dish; or 500µl per 60mm dish**). For cells grown in suspension, spin down cells and lyse directly by the addition of NucleoZOL at **1mL NucleoZOL/10<sup>7</sup> cells**. Lyse cells by pipetting up and down several times. Ensure complete lysis by repeated pipetting.
- 2) **PRECIPITATE PROTEINS/CELL DEBRIS:** Add **200µL** RNase-free water per **500uL NucleoZOL** to the lysate. Shake the sample vigorously for 15 seconds. Centrifuge samples for 5~8min at 12,000 x g at 4°C or room temperature.
- 3) **PRECIPITATE TOTAL RNA:** Carefully pipette RNA containing supernatant from step 2 into a 1.7ml or 2.0ml RNase-free Eppendorf tube(s). Add **1mL isopropanol** per **1mL supernatant (i.e., equal volume)** in order to precipitate RNA. **Make sure tubes are kept on ice**. Centrifuge samples for 5~8min at **12,000 x g** at 4°C or room temperature. **(NOTE: You can stop here by storing the isopropanol-RNA mix at -80°C for days to years).**
- 4) **WASH RNA:** Use **400–600µL 70 % ethanol** when precipitating in 1.5mL or 2mL tubes. For larger tubes, add **500-700uL 70 % ethanol per 1mL supernatant**. Centrifuge the pellets for 1~3 min at 12,000 x g.
- 5) Remove ethanol from the pellet by **carefully pipetting**. Repeat the ethanol washing step. Do not dry the RNA pellet between wash steps.
- 6) **RECONSTITUTE RNA:** Dissolve the RNA pellet in RNase-free water to obtain an RNA concentration of 1~2µg/µL (**i.e., usually dissolve in 20~30µl for one 100mm dish, or 10~15ul for one 60mm dish**). Vortex the sample for 2~5min at room temperature to dissolve total RNA.
- 7) Your RNA samples should be aliquoted and kept at -80°C (**not your own -20°C freezer space**).

**NOTE:** Similar to TRIzol reagent, NucleoZOL reagent contains a very low level of quinidine to inhibit RNase activity. Thus, you should take the following considerations in mind:

- 1). You should try to collect one sample at a time (keep the rest in the incubator). Once you lyse the cells, please collect the lysate and transfer to 2ml-Eppendorf tubes. Please keep the tubes on ice.
- 2). The RNA will severely degrade if the cell lysate is kept in -20°C or -80°C. You should process the samples all the way to isopropanol precipitation/ethanol washes (even though you may leave the isopropanol-filled tubes at -80°C for weeks/months).
- 3). If you plan to store your RNA samples for quite long periods, you can dissolve the RNA samples and perform PC-8 extraction, followed by ethanol precipitation/washes. PC-8 extraction will remove RNase completely.