

Checking Competency of Electrocompetent Bacterial Cells

Adapted from BV's Cookbook, TCH 1/27/02

Titration of electroporation efficiency

| | |
|--------------------------------------|------------|
| CsCl-purified pBlueScript (10 ng/ul) | 1 μ l |
| DH10B electrocompetent cells | 20 μ l |

- 1) Transform cells with Bio-Rad gene pulser at 1.8KV using a 1mm cuvette.
- 2) Add 1 ml of LB at room temperature.
- 3) Shake cells at 37°C for 1 hour in a 15ml tube.
- 4) Make serial dilution of the cells with LB.
10 ng/ml \rightarrow 100 pg/ml \rightarrow 10 pg/ml \rightarrow 1 pg/ml
- 5) Plate 0.1 ml of 1 pg/ml on L-Amp.
- 6) Incubate overnight at 37°C.

NOTE: If 1,000 transformants grow on the plate, efficiency= $1,000 \times 1 \mu\text{g}/0.1 \text{ pg} = 1 \times 10^{10}/\mu\text{g}$ of pBS. Efficiency should be approx. $3 - 10 \times 10^9/\mu\text{g}$ pBlueScript