

NORTHERN BLOTTING GEL TRANSFER AND HYBRIDIZATION

Hongwei Cheng 1/30/02; commented by TCH 1/31/02

I-Buffer preparation:

1. 50mM NaOH/10mM NaCl:
1N NaOH 25ml
5M NaCl 1ml

 Q.S. to 500ml
2. 10x Sodium Phosphate Buffer:
Na₂HPO₄ 160.8gm
NaH₂PO₄ 55.2gm

 Q.S. to 4 L.

II-Pretreatment of gel:

- Wash with ddH₂O x 5-10min.
- Wash with 50mM NaOH/10mM NaCl x 20-30min.
- Wash with 1x Sodium Phosphate Buffer x 10-20min.

III-Gel Transfer:

- cut proper size of filter paper, membrane (Immobilon-Ny, Millipore) and paper towels.
- Set up transfer device (from bottom to top):
glass plate > bridge filter paper (two pieces) > filter paper (two pieces) > gel (face down)
> membrane > filter paper (two pieces) > paper towels > glass plate > weight.
- Transfer in 1x Sodium Phosphate Buffer for 24hrs.
- Remove paper towels and the gel (**you can check on UV box for the presence of residual RNA**). Mark the top of each lane position on the membrane (**NOTE: Always handle the membrane with tweezers or forceps thereafter**).
- Wash the membrane with 1x Sodium Phosphate Buffer, let air dry briefly on a piece of filter paper (**RNA side facing up**).
- Crosslink at energy = 2000uj/cm² for 30 seconds (UVP CL-1000) with RNA side facing up (**NOTE: You'd better record an image of the RNA gel under UV. But you need to have the RNA side facing down this time!**).

IV-Prehybridization, Hybridization and Washes:

- Put membrane into a plastic bag (**using forceps or tweezers**).
- Add 8-10ml pre-warmed QuikHyb (**prewarm at 68°C, & shake well**) into the bag (**make sure the membrane is all wet and well soaked, of course, no bubbles**).
- Seal the bag after expelling the bubbles.
- Prehybridize at 68°C for 20-30min.
- Add labeled probe directly into the prehyb bag. Carefully seal the outmost opening of the bag, then expel bubbles, and seal close to the edge of the blot membrane.
- Mix well by massaging the bag (**also check any leaks**).

- Hybridize at 68°C for 1-2 hrs.
- Wash with 0.3x SSC/0.1% SDS at 65°C for 5-30min (**Note: Monitor the retained radioactivity every 5-10min**).
- Wash twice with 0.1x SSC/0.1% SDS at 65°C briefly (Again, **monitor the retained radioactivity every 5-10min**).

Wrap the membrane, and expose to the x-ray film or a phospho Imager.