

**MOLab Optimized Decalcification Protocol for Decelleularized Bone Matrix (DBM)
(or Bone Disc) Scaffolds**

(Meng Zhang @ 03/12/2020; edited by TCH)

- 1) Fix the retrieved DBM masses in 5-10 volumes of 10% PBS buffered formalin for 12-24h.
- 2) Wash/rinse all fixed samples in slowly running tap water for 30-60 minutes. Avoid rinsing in rapidly running water.
- 3) [**Optional:** You can do microCT imaging at this stage].
- 4) Decalcify samples in ~20 volumes of **10% nitric acid** for **5~7h at 37 °C** water bath with gentle agitations [**Note:** you should check the decal status every 30min after 5h to avoid excessive decalcification in nitric acid].
- 5) [**Optional:** Rinse samples in slowly running tap water for 30-60min; leave samples in our regular Cal-Ex II decal-fix solution for another 12-24h].
- 6) Once the decalcification is complete, rinse samples in running tap water thoroughly overnight, followed by paraffin embedding using the automatic tissue processor in J-337.

NOTES:

- 1) Samples should not be clustered together nor sit at the bottom of container in order to provide for complete decalcification;
- 2) To check the end-point of decalcification, simply insert a pin/syringe needle to see if they are easy to be pierced.