MOLab Standard Protocol for Preparing RT-PCR cDNA Stocks for qPCR

(Lan Zhou 11/02/00, Yang Bi 02/24/08, updated by Yi Shu @ 03/30/2025; commented by TCH)

NOTE:

- A) Yang has recently compared the RT enzyme from NEB (Cat# M0253S) vs. Invitrogen's Superscript II, and found no significant difference between those enzymes. Thus, NEB RT enzyme is the recommended choice in our lab. However, we do not use any manufacturer's protocol!
- B) After numerous pre-testing runs by former/current lab members, we decide that the RT enzyme in each reaction can be reduced to as low as **0.25µl** (the original suggested volume was 1.25µl).
- C) These conditions work well for most regular RT-PCR cDNAs prepared for real-time PCR assays, especially if your RNA is isolated from cultured cells and/or fresh tissue samples. However, if you are concerned about the yield and quality of your RNA samples, you are always recommended to check a couple microliters of your RNA samples on RNA gel (see Protocol C12).
- D) Make sure you use enough total RNA for each RT reaction. Typically, a subconfluent T-25 flask yields about 5-10µg total RNA, and a subconfluent 100mm dish yields about 10-20µg total RNA. Please keep in mind that the RNA yield will be much lower if the cells are cultured in low serum medium 8-12h prior to RNA isolation.
- E) In our TRIzol or Nucleozol RNA isolation protocol, we do not routinely use RNase-free DNase to digest away potential genomic DNA. Thus **-RT tube** is set up in some cases.

RT-PCR PROTOCOL

- 1. Turn on two heating blocks: one at 70° C and one at 42° C.
- 2. Make Hexamer mix-RNA [prepare half of the volume if -RT tube is not set up]:

| Total volume | 12 µl | |
|---------------------|-------|--|
| Total RNA | 10 µl | |
| Hexamer (0.5 µg/µL) | 2 µl | |

- 3. Incubate @ 70°C for 5 minutes.
- 4. Prepare RT mix [prepare half of the volume if -RT tube is not set up]:

| Total volume | 8 2 ul <mark>(~8ul /RT F</mark> | |
|--------------------|---------------------------------|--|
| ddH2O | 3.5 µl | |
| M-MuLV | 0.2µl | |
| 5mM dNTPs | 2 µl | |
| 10 X M-MuLV buffer | 2.5 µl | |

5. Prepare RT Reactions:

| | +RT tube | |
|-----------------|----------|--|
| Hexamer mix-RNA | 12.0 µl | |
| RT mix | 8.0 µl | |
| Total volume | 20µl | |

7. The cycling program of RT-PCR:

42°C X 60 minutes 95°C X 1 minute (for killing RT; RT interferes with Taq DNA Pol) 4°C Hold or kept at -80°C freezers

- Add 80µl ddH₂O to get 100µl original cDNA, and then take 10µl out, add 40µl ddH₂O to dilute into 5 times (total volume = 100µl, aka., RT-PCR cDNA stock), and then use the diluted cDNA for real-time PCR. *It's a good practice to aliquot the initial 100µl into multiple aliquots. Keep the aliquots at -80°C.*
- 9. For most real-time PCR reactions, the cDNA mix needed to be further diluted 5 to 100 times, depending on transcript abundance of the gene of your interest. *The ideal Ct or Cq values for internal reference genes should be between 15~20 cycles.*